



LC880 Experiment Controller User Guide

Documentation for the LC880™ Controller
and Trigger™ Control Software v. 5.0

Trigger software ©1996--2008 LabSmith.

This manual ©2008 LabSmith. No part of this document may be reproduced or distributed without the consent of LabSmith.

TABLE OF CONTENTS

INTRODUCTION	5
Parts List	5
Specifications	5
Basic Functions	6
How You Will Use the LC880	6
Installing Trigger Software	7
Setting Up the LC880	7
Field Upgradeable Software	8
Getting Help	8
OVERVIEW OF LC880 CONTROLS	9
Back Panel Inputs and Outputs	9
Front Panel Controls	10
TRIGGER SOFTWARE INTERFACE	11
Main Screen	11
General Settings	13
Disable/Enable I/O	13
Lock/Unlock Panel	14
Time base settings	14
CHANNEL PROPERTIES	15
Channel Properties Screen	15
Quiescent States	16
Logic	17
Logic Assignments	18

Parsing	19
Sending Logic	19
Experimenting With Channel Logic	20
Rules for Logical Assignments	21
Predefined Variables	21
Global Variables	23
Timing Modes	24
“Apply” versus “Program All Channels”	25
Fixed Output Mode	26
Passive Mode (Output = Input)	27
Clocked Pulse Stream Mode	28
Delayed Pulse After Trigger Mode	30
Important Delay Timing Considerations	35
Counter/Timer Mode	39
Dynamically-Delayed Pulse Mode	41
Toggled Output Mode	46
Dynamic Delay Compensation Mode	47
RUNNING EXPERIMENTS	50
Storing and Recalling Experiment Files	50
Front Panel Controls and Indicators	50
Channel LEDs	50
ALL Channel LED	52
Selecting a Channel from the Front Panel	52
Performing a Front Panel Function	53
Enable/Disable ALL	53
Enable/Disable Channel	53
Trigger a Channel	54
Reset	54
Store/Recall Settings Profiles	54
GLOSSARY	56

INTRODUCTION

The LC880 is a fully programmable logic and timing controller for coordinating and synchronizing lab equipment and physical experiments. LC880 provides clocks, counters, triggers, and many other useful functions, with 10 ns resolution and 100 ps accuracy.

As you are well aware, controlling and synchronizing experimental equipment are often the most tedious tasks you'll face in the lab. Experimenters regularly build their own controls from scratch or piece together systems from the limited offerings on the market. The LC880 simplifies your control setup, so you can go about the business of gathering data, and making knowledge.

The LC880 consists of two components: Trigger software and the LC880 experiment controller. Trigger lets you design and run experiments from a simple graphical interface. The LC880 controller provides the connections and control for your experimental devices.

Parts List

Your LC880 package should include the following items:

1. the LC880 controller
2. AC power cable
3. Trigger software installation CD
4. 9-pin RS232 cable
5. User's manual.

If any parts are missing or damaged please contact LabSmith support.

Specifications

Trigger software will run on any PC-compatible computer running Microsoft® Windows® 2000 or XP or later. Drivers for National Instruments® LabView™ software are also available for download. See <http://www.labsmith.com> for more information.

The LC880 includes eight TTL-compatible inputs and eight TTL-compatible outputs. Voltages in the range 2.5–5 V are TTL “highs” and voltages in the range 0.0–0.2 V are TTL “lows.” Intermediate voltages should not be present on inputs and outputs except during transitions. The inputs and outputs are protected against connection to 115 VAC. When disconnected the inputs are read as TTL “highs.”

The LC880 can be connected to/ disconnected from the PC at any time without damage to either unit. This ability allows you to program to LC880 at your computer, then to take the unit to the experiment area for stand-alone use.

Basic Functions

The LC880 provides eight functional **channels**, or internal signal paths.

At any time, each channel can be assigned one of eight **timing modes**. Timing modes include clocks, counters, delayed triggers, etc. You can set various **properties** for each timing mode—for example, the delay length on a Delayed Trigger. All modes are aligned to a precise, 100 MHz internal **clock or an external clock**, allowing resolution down to 10 ns.

Each channel also includes **logic** settings which specify the conditions that will trigger the channel. A channel can be triggered by any logical combination of signals from the eight inputs and/or the outputs of other channels.

The timing properties and logic that you set up in Trigger can be transferred and stored in the memory of the LC880. You can store up to 64 variations of these properties in **settings profiles** in the LC880. You can then switch between settings profiles to switch experiments, compare different timing scenarios, or simulate different real world conditions.

All settings and setting profiles for a particular experiment can be stored as an **experiment file** on your computer. Experiment files bear a .trg extension.

How You Will Use the LC880

In most cases, you'll use the LC880 as follows:


1. Using the Trigger software, open an existing experiment file, or create a new file.
2. In Trigger (or LabView), program the properties and logic for each channel and for the experiment as a whole.
3. Send the properties and logic to the LC880's memory.
4. Test your program by simulating triggers and examining the outputs using the LC880's channel LEDs and/or an oscilloscope or other diagnostic tool.
5. Store variations of channel properties in one of the LC880's 64 settings profiles.
6. Plug your experimental equipment to the LC880.
7. Run the experiment, using the LC880 by itself or leaving it tethered to the PC. Switch between settings profiles to explore various timing conditions.

Installing Trigger Software

Insert the Trigger Software disk in your CD-ROM drive and execute the setup.exe file.

Setting Up the LC880

1. Plug the power cable from the back of the LC880 into a 100–250 VAC power outlet.
2. Connect the 9-pin cable to the back of the LC880 and to a free serial port on the PC.
3. Press the Power button on the front of the LC880.
4. From your computer, run Trigger.exe to launch the Trigger software.
5. From within Trigger, select View>Options>Com Port. The Serial Communications screen will open.



6. Select the port to which the LC880 is connected. The LC880 will automatically detect the baud rate, which is typically set to 38400. If the RS232 cable running to your LC880 is greater than 30 m long, you may need to lower the baud rate to achieve reliable communication.

Field Upgradeable Software

The firmware in the LC880 is field upgradeable, allowing access to new timing modes and capabilities as they become available.

To upgrade the firmware, download the latest version from www.labsmith.com. Then choose “Update Firmware” on the Program menu to load data from that file.

Getting Help

This guide is your main source for information on operating the LC880 and Trigger software. The guide is also available in an Adobe® Acrobat® Portable Document Format (pdf) file for electronic viewing.

While running Trigger, press F1 to view this document on-line. Or, press SHIFT-F1, then select a button or menu option to view help topics for a particular feature.

The LabSmith web site (www.labsmith.com) also includes helpful applications notes and other technical information to help you use the LC880.

If you are unable to find the information you need, please contact our technical support team at (925) 292-5161, or e-mail us at info@labsmith.com.

If you need to contact support, please let us know your LC880 serial number (located on the bottom of the unit) and the version of software you are using. **To determine the software version number**, select “About Trigger” from the Help menu.



OVERVIEW OF LC880 CONTROLS

The LC880 is the control unit and interface for your experimental equipment. You will program it using the Trigger software. Once you've designed, tested and downloaded your program, you can control the LC880 from the computer, or you can run it in stand-alone mode.

Indicators and controls on the LC880 front panel will help you run your experiments. You'll learn all about these functions in the *Experimenting* section later. For now, we'll give you a short tour of the inputs, outputs and controls so you can work with the unit while you learn to use the Trigger software.

Back Panel Inputs and Outputs



The back panel includes the following connectors:

- an internally-fused jack for a standard 100–250 VAC power cord
- a 9-pin RS232 jack to connect to the Com port of your PC
- eight BNC channel input jacks (in1–in8), normally high (5 V)
- eight BNC channel output jacks (outA–outH), normally low (0 V)
- a BNC input jack for an external clock source.

The inputs and outputs will accept any device into that sends or receives TTL-compatible signals.

Front Panel Controls



The LC880 front panel controls begin with the **Power button** and **Power LED** in the upper left corner. When you first press the Power button, the LED will flash green and red while the system runs its initial diagnostics. It will then glow steadily: if it glows a steady green, the front panel is unlocked and accepting button pushes; if it glows red, the front panel is locked (see *Trigger Software Interface: General Settings* below). Even the power switch will be deactivated, safeguarding against mis-setting or tampering.

The **channel LEDs** offer useful information about the state of each channel. A green light typically indicates that the channel is ready to receive a trigger. A red light typically indicates that the channel is in the process of triggering. When a channel indicator is off, the channel is disabled (see *Running Experiments: Enable/Disable Channels*). More will be discussed in the *Timing modes* section.

The **All Channel LED** shows the status for all channels. A green LED indicates that the LC880 will respond to external inputs and that the external outputs are active. A red LED indicates that the LC880 will not respond to external inputs and that outputs are in a fixed and frozen logical state.


The **Channel and Function Buttons** work as pairs—you will select one of each to perform any function. In the base state, all buttons glow dimly, indicating that they're ready for your instruction. Once you press a channel or function button, only the available options will be illuminated.

TRIGGER SOFTWARE INTERFACE

Trigger has been designed exclusively for setup and control of the LC880. You'll use Trigger to set the timing properties and logic for each of the eight channels and to set up overall properties for the experiment.

Main Screen

When you run Trigger, you'll be greeted by the Main screen:



NOTE: At start-up, the Main screen will show the default program, which sets all channels to “delayed trigger” mode.

From this screen you'll control or access all of Trigger's functions. Most Trigger functions can be accessed in two ways by selecting an option from a pull-down menu or by pressing a button on the Main screen.

1. Menu options:

The **File menu** lets you open, close and save experiment files (.trg files).

The **Edit Menu** lets you cut, copy and paste text as you edit your program.

The **View menu** lets you set up the Trigger interface. Choose to display or hide the tool bars and status bar, select the font used in Trigger, and set up communications with the LC880.

The **Channel menu** lets you open the screens for “General Settings” and for each of the eight channels.

The **Program menu** includes options for communicating with the LC880. You can send logic, send properties, reset all channels, or enable/ disable the LC880’s inputs and outputs.

2. The buttons are as follows:



The **General Settings** button opens the **General Settings** screen.

The **Channel** buttons open the **Settings** screens for all channels.



The **Program** button updates the LC880 with all current settings.




The **Reset** button resets all channels of the LC880. Retriggering Counters and delay settings will be reset as well.

The **Main** window panel shows a text summary of all of the settings in the current experiment. This summary may be printed for review and documentation purposes. The text is not freely editable and must be modified through the dialogs. Double-clicking on selected parts of the summary opens the dialog needed to edit the text.

At any time, the LC880 control bar can be hidden or displayed by clicking View>LC880.

General Settings

The “General Settings” screen lets you control functions that affect all channels of the LC880. To reach this screen, click the “General Settings” button, or choose “General” from the Channel menu.



Disable/Enable I/O

On occasion you will need to abort an experiment or quickly place all outputs to your equipment in a “safe” state. Each of the tabs of the General Settings screen includes buttons that let you “disable” and “enable” all inputs and outputs to the LC880. With I/O disabled, input signals are ignored. Each output will be held at either a logic low or logic high, depending upon its channel properties. We’ll refer to the “safe” voltage that a channel produces when disabled as the “quiescent state” (see *Channel Properties: Quiescent States*). To restore channels functions, click the “Enable I/O” button. These functions are nonvolatile: the enabled and disabled states are restored on power cycle.

You can also disable or enable inputs and outputs from the LC880 front panel (see *Experimenting: Disable/Enable All*).

Lock/Unlock Panel

Lock or unlock the front panel buttons using the “Lock Panel” and “Unlock Panel” buttons. This functionality is handy, for example, in an embedded control application where the LC880 should function only in stand-alone or computer-controlled modes. When the panel is locked, the LED next to the power button glows a steady red and all buttons are dark. When the panel is unlocked, the power LED glows green and the front-panel buttons are illuminated.

Like the enable/disable I/O functions, these settings are nonvolatile and are restored on power up. Because the power button is disabled, the LC880 must be disconnected from AC power to turn it off, if the front panel is locked.

Time base settings

Each pair of channels of the LC880 can independently utilize one of four different time bases so you can control and measure events which occur a few nanoseconds apart, or a few hundred seconds apart.

In “**100**” MHz time base mode (the default), the LC880 aligns signals to an internal 100 MHz clock, for timing resolution of 10 ns. Use this “high resolution” mode when the events of your experiment occur in rapid succession (milliseconds or less).

“**100/64**” MHz time base mode allows the LC880 to control events that occur over longer time periods (milliseconds to minutes). This is called “low resolution” mode in the following documentation.

In “**Ext.**” time base mode, the channels use an external clock input for the time base. This input must be TTL-level between 1 and 100 MHz applied to the External Clock (ext clk) input on the back of the LC880. This function allows the use of high-precision or variable-frequency time bases. It also allows multiple LC880s to be synchronized to the same base.



NOTE: The “External” time base cannot be produced by a channel that uses the same time base as an input.


“**Ext./64**” time base mode divides the frequency of the signal applied to the external clock input before using it as a time base. If the external clock is being used as a time base, the frequency entered into the edit box is used to calculate the timing settings. If this setting is 10% higher than actual, the channels using this time base will be 10% long on delays and duration.

CHANNEL PROPERTIES

Channel Properties Screen

Each channel has a Settings screen in which you'll select the timing mode, set the timing properties, and establish the logic that will determine when the channel will trigger.

To open the Settings screen for a channel, select the channel from the Channel menu, or click a channel button (A–H) on the Main screen. You can open more than one channel screen at a time.



Channel properties are organized by tabs. The tabs that are relevant for a particular timing mode appear automatically when that timing mode is selected. For convenience, all settings are retained and saved, including those for timing modes that are not currently selected.

Four additional buttons appear next to the tabs:

- **Disable** temporarily disables the channel. The channel remains fixed in its “quiescent state” (see *Quiescent States* below). The front panel LED for that channel will be unlit (See *Running Experiments: Channel LEDs*).
- **Enable** restores channel functions following a temporary disable. The front panel LED for that channel will glow to indicate the channel's status.
- **Reset** reinitializes the internal trigger counters and sweep timing for the current channel. This function can be called simultaneously for all channels by pressing the “Reset” button on the main tool bar.
- **Force Trigger** lets you simulate a valid input signal for the channel (see *Channel Properties: Logic* for more). This button, which is only available with functions that *can* be triggered, is useful both for troubleshooting your channel settings and implementing software control of experiments. You can also force-trigger a channel from the LC880 front panel (See *Running Experiments: Force Trigger*).

Quiescent States


When you disable a channel, its input will not change the output. The “quiescent state” of a channel is the setting of its output when the channel is disabled. You can set each channel’s quiescent state to be a TTL high (3.5–5 V) or low (0–0.2 V) to provide a safe signal to the output device if you need to abort an experiment or stops its progress.

The method for setting the quiescent state depends on the timing mode. Most modes include an “Invert Output” option which sets the output normally High.

Logic

A channel can be triggered by any combination of signals from the eight inputs and/or other channel outputs. Channel **logic** determines the combination of signals that triggers the channel.

The diagram below shows the architecture of the LC880. Note that any inputs and outputs can trigger a channel without external cabling.



To define when a channel will trigger, open the channel's Settings screen then click the "Logic" tab. This tab only appears for timing modes which respond to logic. You will enter a logic assignment for the channel in the large white box.



NOTE: Some timing modes do not respond to trigger logic. For example, a channel set to operate as a clock will continue to run regardless of what happens to its input.

Logic Assignments

The logic processor assigns the combination of external inputs (in1–in8) and channel outputs (outA–outG) that will trigger each channel. Logic assignments always appear in the following format:

[the channel to be triggered] = [the conditions that will trigger it];

External inputs and channel outputs feed into the logic processors; therefore, they appear only on the right side of assignments. Logic processor outputs feed the timing channel inputs inA–inH; these terms appear only on the left side of assignments.

Here is an example:

inB = in2;

This assignment means, “Channel B will be TRUE when Input 2 goes TRUE (reaches ~ 5 V).”

The assignment:

inB = not in2;

means “Channel B will trigger when Input 2 goes FALSE (drops below 0.2 V).”

You can employ compound arguments in logical assignments, such as:

inB = not outC and in2;

This means, “Channel B will trigger when the output of Channel C is FALSE *and* Input 2 is TRUE.”

This example is more complex still:

inB = (in1 or in2) and (outC or outD);

Here, Channel B requires a TRUE signal from Inputs 1 or 2, combined with a TRUE signal from Channels C or D, in order to trigger.



NOTE: LC880 inputs are TRUE by default if nothing is connected (i.e., the input is normally weakly pulled high). Grounding an input switches its state to a logical FALSE.

It is important to distinguish between the eight physical input jacks, in1–in8, and the eight timing channel inputs. Any combination of the *physical* inputs can be used to trigger the *channel* inputs. The following assignment is incorrect:

```
outA = in1;
```

because it attempts to trigger a *channel* output rather than a channel input, inA–inH

This assignment is also invalid:

```
in1 = in2;
```

because it attempts to assign a value to External Input 1.

Finally, this assignment is invalid:

```
inB = inA;
```

because it attempts to assign a value from one channel input to another.

Parsing

When you have completed a logical assignment, click the “Apply” button. Trigger will check your syntax and logic and will report any mistakes to you by line number.

Sending Logic

If you change channel logic you must transfer those changes to the LC880 to implement them.



To transfer logic to the LC880, click the “Set All Channels” button , or choose Program>Program All Channels.

While the logic is transferring, the Power LED on the LC880 will turn red. This transfer typically takes only a second, though it may take slightly longer for complex programs.

After the transfer, the logic settings reside in the LC880, even if you turn the unit off or disconnect it from the PC.

Experimenting With Channel Logic

This brief example will show you more about channel logic:

1. Connect a normally-open switch to Input 1 and another to Input 2.
2. Connect an LED to Output A.
3. Open the Channel Settings screen for Channel A.
4. Choose “Output = Input” from the “Timing Mode” list. In this mode, the output of Channel A is always the same as its input.
5. Select the “Passive” tab and choose “Output = Input.”
6. Select the “Logic” tab and enter:

$\text{inA} = \text{in1};$

7. Click the “Set All Channels” button.

The LED at Output A should be glowing, indicating that the output is TRUE, or HIGH (since LC880 inputs are normally HIGH). When you close Switch 1, the LED will turn off, indicating that the output is FALSE, or LOW.



NOTE: The Channel A LED on the front of the LC880 glows green while the channel output is LOW, or red when the output is HIGH. But don't be confused! The channel LED does *not* always indicate the channel's output state—only in certain Timing Modes. You'll learn more later in the *Timing Modes* section.

8. Next, change the logic assignment to:

$\text{inA} = \text{not in1};$

9. Click the “Set All Channels” button. The LED at Output A will now remain off until you close the switch.

10. Change the logic assignment to:

$\text{inA} = \text{not in2};$

11. Click the “Set All Channels” button. The LED at Output A will now remain off until you close Switch 2.

12. Lastly, change the logic assignment to:

```
inA = not in1 and not in2;
```

13. Click the “Set All Channels” button. The LED at Output A will remain off until both switches are closed.

Rules for Logical Assignments

The following rules apply to logical assignments:


- Logic is case-sensitive. Use lowercase and capital letters as in the samples above.
- Spaces matter. Use spaces per the samples above.
- If a logical assignment consists of more than one line, each line should end in a carriage return.
- The last line of a assignment must end in a semicolon.
- You can add comments to annotate your logic assignments. Single line comments begin with //. Multi-line comments begin with /* and end with */.
- Many symbols can be written using C-Style programming commands. For example, “not in1” can be written “!in1”.

The following symbols can be used in logical assignments:

<u>Symbol</u>	<u>Shorthand</u>
not	!
and	&
or	(vertical line)
xor	#
HIGH	<i>none</i>
LOW	<i>none</i>
()	<i>none</i>
;	<i>none</i>

Predefined Variables

Trigger includes a number of Predefined variables to speed logic coding. To view the variable names that are recognized by Trigger, click “Predefines” on the “Logic” tab.



- inA–inH: the input of channels A through H. These signals only appear internally at the input of the timing channel; they are not connected to the rear BNCs. These symbols can only be on the left-hand side of an assignment
- in1–in8: inputs 1 through 8 from the rear BNC connectors.
- outA–outG: outputs of channels A through G that appear at the rear BNC connectors.



NOTE: The output of channel H is not internally connected to the logic processor. It cannot be used in logic assignments unless it is externally connected to another channel input.

- LOW, FALSE: trigger input always low.
- HIGH, TRUE: trigger input always high.

Double click a variable to copy it to your Logic assignment. You can also select the variable then click “OK.”

Global Variables

“Local Variables” are those that are used directly in logic assignments, such as “inA” or “out1.” Logical assignments that are used repeatedly in your experiment can be stored as “global variables.” With global variables, complex logic is easier to read and debug.

Several channels in your experiment may need to be triggered by the same event. For example, Channels F, G and H must trigger when:

(in1 or in2) and (outC or outD).


Rather than typing that logic each time, you may define it as a variable:

event1 = (in1 or in2) and (outC or outD);

Now you can use “event1” in the logic for Channels F, G and H:

inF = event1;

To view the global variables, click the “Globals” button on the “Logic” tab. The “Globals” screen will open:



To create a new global variable, type the variable's name followed by its definition. In this example we've defined variables named "arm" and "fire." Be sure to end your definition with a semicolon. Click "Apply" to save the variable.

Timing Modes

In the previous section you learned that logic defines how channels are triggered. In this section you'll use the LC880 timing modes to determine how the channels will respond.

Each channel can be configured to operate in one of eight modes:

- Fixed output
- Output=Input (Passive Mode)
- Delayed Pulse After Trigger
- Clocked Pulse Stream
- Validated Trigger Input (Noise Suppression)
- Dynamically Delayed Pulse
- Counter/Timer
- Toggled Output.

When a timing mode is selected, the relevant property tabs will be displayed.



NOTE: All timing property tabs include a space for Notes. These notes are saved and appear in the Main Screen program summary area. It is good practice to document your settings extensively. As with all programming, good documentation will simplify the process of creating, debugging and running your experiments.

For each timing mode you'll be able to set up a simple example to view the output of the channel.

To set up for the experiments later in this section:

1. Make sure that the LC880 is connected to the PC.
2. Plug a normally open push-button trigger switch into Input 1. The button will represent our triggering "event."
3. Open Channel A's Setting screen, and select the "Delayed pulse after trigger" timing mode (you'll learn about Delayed Trigger mode shortly).

4. Select the “Logic” tab. Enter the logical assignment:

inA = not in1;

Then click OK. With this logic, Channel A will trigger when you depress the trigger button.



5. Now transfer the logic to the LC880 by clicking the “Set All Channels” button
6. Plug an oscilloscope to Output A to observe the channel output. Remember, the output voltage will vary from 0 to 5 V.

“Apply” versus “Program All Channels”

Each properties screen includes an “Apply” button which, when clicked, will send the current timing properties to the channel on the LC880.



NOTE: If you change a channel's properties, closing the Settings screen will have the same effect as clicking “Apply.” Clicking the “OK” button will also have the same effect.


Logic settings are sent and saved in the LC880, but any changes are not implemented on the

LC880 until you click “Set All Channels” button  on the Main screen or choose Program>Program All Channels.

“Program All Channels” is the same as clicking “Apply” for all channels and then instructing the LC880 to update its logic programming.

The “Cancel” button cancels any edits made since the last time the “Apply” button was clicked and closes the dialog window.

Fixed Output Mode




In this mode, the channel output remains in the selected state (Low or High), regardless of the input.

The front panel LED for the channel will glow green when the output is low or red when the output is high.

To try fixed output mode:

On the "Properties" tab of Channel A, choose "Fixed Output" from the pop-up. Choose "Low" or "High" as the output state, then click "Apply" to save the settings. At Output A, the signal will remain a steady TTL logic low or high, depending on your selection.

Passive Mode (Output = Input)




In passive mode, the output signal equals the input signal or its logical inverse.

The front panel LED for the channel will glow green when the output is Low or red when the output is high.

To try passive mode: On Channel A's Settings screen, select the "Output=Input" timing mode. In the Function box, choose "Output=Input," then Click "Apply" to save the settings.

When the trigger button is not depressed (i.e., the input signal is low), the signal at Output A will be a TTL low. Now, depress the trigger button at Input 1—the signal at Output A will go high. If you select "Output=!Input (invert), the channel output will remain high until you depress the trigger.

Clocked Pulse Stream Mode



In clock mode, the output of the channel is a steady stream of pulses. You set the pulse duration and delay between pulses. A clock channel's input does not affect its output.

The front panel LED for a clock channel will glow green when the output is low or red when the output is high.



NOTE: For a high-frequency clock the LED will alternate between red and green so quickly that it will appear as a continuous pale green, yellow or orange light.

As mentioned earlier, the LC880 can operate in two internal clock modes: “100” MHz high-resolution and “100/64” MHz low-resolution (see *General Settings: Time Base Settings*). In high resolution mode, the clock pulse and duration can range from 0.48 μ s to 52.4 ms. In low resolution mode, pulse and duration can range from 30.1 μ s to 3.35 s.

To set the duration of the high state, select the timing units (seconds, microseconds, milliseconds or nanoseconds) from the pop-up, then enter a duration in the “Clock high duration” box.

To set the duration of the low state, select the units, then enter the duration in the “Clock low duration” box.

A clock’s quiescent state is normally low (See *Channel Properties: Quiescent States*). Check the “Inverted output” box if the quiescent state for the channel should be high instead.

To try Clock mode:

1. On the General Settings screen, select “100/64” (low resolution mode) for Channels A and B.
2. On Channel A’s Properties tab, choose “Clock” as the mode, select an “on” duration of 100 ms and an “off” duration of 500 ms. Click “Apply” to save the settings.


The signal at Output A will go high every 500 ms and will remain high for 100 ms each time.

To test faster clock settings, select “100” (high resolution) mode under General Settings, then set duration and delay times of less than 52 ms.

Delayed Pulse After Trigger Mode

Upon receiving an input signal, a delayed trigger waits a specified delay time, then sends a pulse to the channel output. This mode has several property tabs.

Delay



On the Delay tab you'll specify:

- Delay after trigger** The delay between a valid trigger and an output signal
- Duration of pulse** The duration of the output signal
- Until Reset** The output will remain high after a valid trigger, until the channel is reset


Inverted Output Check the box to set the output high until a valid trigger.

The range of settings you can use with the delayed trigger depend upon the current timing mode:

	100 (high res mode)	100/64 (low res mode)
Delay range	20 ns to 21.47 s	1.28 μ s to 1374 s
Duration range	7.7 μ s to 21.47 s.	493 μ s to 1374 s
Delay / Duration resolution	10 ns	640 ns


Triggering

This tab lets you control how the channel will trigger:



Normal	The channel can retrigger an unlimited number of times
Skip N triggers	The channel triggers only after the specified number of input pulses (between 0 and 2,147,483,647).
Retrigger N times	Limits the number of retriggers allowed between resets (between 0 and 2,147,483,647). To allow a channel to be triggered again after reaching the limit, press the “Reset” button on the LC880 or choose Program>Reset Channels in Trigger.
Trigger on rising input	Sets the channel to be edge triggered. In other words, an input signal that stays high triggers the channel only once.
Trigger on input high	Allows the channel to be repeatedly triggered while an input is high. This setting could be used for generating a low-precision gated clock, for example. However, the latency time between retriggers is indeterminate (typically around 10 μ s) so a clock made this way is not suitable for applications requiring high precision.
Alert host computer trigger	Sends a “ping” to the PC upon completion of the trigger following cycle.

Sweep



To facilitate time-delay scanning experiments, the delay of the pulse relative to the trigger can be “swept” or incremented by a specified value each time the channel is triggered. The sweep value must be positive. The increment resolution is 5 ns. However, the actual delay is rounded to the nearest 10 ns at full internal clock. For example if the increment is set to 5 ns, the actual delay is incremented by 10 ns every second trigger event.

The front panel LED for the delayed trigger will glow green until an input signal is detected. The LED will be unlit for the delay time and will glow red for the duration time. What happens next depends upon the channel settings:

Retriggering enabled?	Edge triggered?	The channel LED...
no	yes or no	remains red (no more triggers accepted)
yes	no	glows green (waiting for next trigger)
yes	yes	is off; glows green once input goes low.

To try the delayed pulse mode:


Under General Properties, choose 100/64 (low resolution) clock mode for Channels A and B.

1. In Channel A's properties screen, select "Delayed Trigger."
2. On the Delay tab, enter a delay of 1 s and duration of 500 ms.
3. On the Triggering tab choose "Retrigger N Times" and set N=0. Choose "Input high (level trigger)," then click "Apply" to save the settings.
4. Click the trigger pushbutton. One second after you push it, the output of Channel A will go high for 500 ms.
5. The channel LED remains red, indicating it is not ready to trigger. Push the button again. You'll see that you cannot retrigger the channel. Click the "Reset" button on the LC880—the LED glows green and you can once again trigger the channel.
6. On the Triggering tab, choose "Retrigger N times," set N=5, then click "Apply." Depress and hold the trigger button. As long as the trigger is depressed (i.e., the input signal is high) the channel will continue to trigger until the limit is reached.
7. Again on the Triggering tab, select "Rising Input (Edge Trigger)" then click "Apply." Depress and hold the trigger button. The channel will trigger only once, no matter how long you hold the button. Release the button and depress it again, and the channel will be triggered once more (switch bounce may produce triggers when the button is released).
8. On the Delay tab, check "Inverted output," then click "Apply." The output will be the inverse of that described in Step 4.
9. On the Sweep tab, set the sweep to 0.5 s per trigger, and on the Delay tab select Normal. Click "Apply." The time delay before the output pulse will increase with each trigger.

Important Delay Timing Considerations

Once the LC880 has begun a delay sequence, it ignores incoming trigger pulses. There is also a delay of approximately 10 μs after the completion of a delay sequence while the channel resets for the next trigger. This effect can be used to skip triggers, but if misunderstood it can cause unexpected results.

For example, in the diagram below, two trigger pulses arrive 100 μs apart. The channel is functioning in Delayed Pulse Mode, with a delay of 50 μs and a duration of 45 μs pulse. It would: receive the first pulse; wait 50 μs ; send a 45 μs output pulse; and trigger again at the second pulse. But the extra 10 μs delay at the end of the sequence will cause the channel to ignore the second trigger pulse.



Validated Trigger Input Mode (Noise Suppression)

The screenshot shows a dialog box titled "Channel B" with a close button (X) in the top right corner. The dialog is divided into three tabs: "Logic", "Validate", and "Triggering". The "Validate" tab is currently selected. Inside the "Validate" tab, there are several controls:

- A "Timing mode" dropdown menu set to "Validated trigger input".
- Three sub-tabs: "Logic", "Validate", and "Triggering".
- A "Minimum duration of valid trigger (delay)" section with a text input field containing "30" and a unit dropdown menu set to "us".
- A "Duration of pulse" section with a radio button selected for "30" (with a text input field containing "30" and a unit dropdown menu set to "us") and another radio button for "Until reset".
- An unchecked checkbox labeled "Inverted output (pulse from high to low)".
- A "Notes" section with an empty text area.

On the right side of the dialog, there are four buttons: "Enable", "Disable", "Reset", and "Force Trigger". At the bottom of the dialog, there are three buttons: "OK", "Apply", and "Cancel".

The noise suppression function eliminates false triggers due to noise on the input signal by requiring the signal to remain high (or low) for a specified duration. If the signal drops out during that period, the output will remain low and the channel will look for the next input signal. Once an input signal is deemed “valid,” a pulse of user-specified duration will be sent to the channel output.

Unlike an analog filter which “smears” the signal edge, this noise suppression function digitally blocks false triggers, retaining the sharp edges on the real trigger signal.

As with Delayed Pulse After Trigger mode, you have additional controls over the channel.

On the **Validate** tab, choose:

- Minimum duration** The minimum length of an input pulse that will trigger the channel
- Duration of pulse** The duration of the output signal
- Until Reset** The output will remain high after a valid trigger, until the channel is reset
- Inverted Output** Check the box to set the output high until a valid trigger.

On the **Triggering** tab choose:

- Normal** The channel can retrigger an unlimited number of times
- Skip N triggers** The channel triggers only after the specified number of valid input pulses (between 0 and 2,147,483,647).
- Retrigger N times** Limits the number of retriggers allowed between resets (between 0 and 2,147,483,647). To allow a channel to be triggered again after reaching the limit, press the “Reset” button on the LC880 or choose Program>Reset Channels in Trigger.
- Trigger on input high** While the input is high, channel can be repeatedly triggered.
- Trigger on rising input** Sets the channel to be edge triggered (i.e., an input signal that stays high will only trigger the channel once).


The front panel LED for the noise suppression channel will glow green until an input signal is detected. It will then glow red until the duration period has elapsed. At this point, the LED will respond as follows:

Retriggering enabled?	Edge trigger?	LED
no	yes or no	remains red (no more triggers accepted)
yes	no	glows green (waiting for next trigger)
yes	yes	is off; glows green when trigger goes low

To try validated trigger mode:

1. On the General Settings screen, select 100/64 (low-resolution) mode for Channels A and B.
2. In Channel A's properties, select "Validated Trigger Input" mode.
3. On the Validate tab, enter a "Minimum duration" and "Duration of output pulse" of 1 second each.
4. On the Triggering tab, select "Retrigger N times," and set N=50. Click "Apply" to save the settings.
5. Now, depress the trigger button. If you hold the button down for more than one second (the minimum duration), the output of Channel A will go high. If you release the button before one second elapses, the output will remain low, and the channel will wait for another input pulse.

Counter/Timer Mode



The counter/timer mode determines the time delay between two separate input signals and reports it back to the PC. You can view the counter output in the Properties tab.

Use the channel's logic to select the input signals that will be compared by the counter (see *Channel Properties: Logic* above). For example, if the logic reads "in1 and in2;" the counter will determine the time delay between a pulse at Input 1 and the next pulse from Input 2.

Normally, the counter will report a value each time a pair of pulses is measured. To retain the current measurement, check the "Hold measurement" button; the counter will not update the delay display until this button is cleared.



NOTE: To switch from counter mode, first click the Disable button and ensure the channel LED is unlit, indicating that the channel is disabled. This will help you avoid unexpected behaviors when the counter is transmitting measurements rapidly (e.g., >1000 measurements per second).



NOTE: When you first run the LC880, channels that are configured as counters are initially disabled and must be enabled either by clicking the “Enable” button on the channel’s Properties screen, or by using the front panel (See *Experimenting: Disable/Enable Channels* below).

The front panel LED for the counter will be green until the first trigger, red until the second trigger and off while the count data is reported and the channel waits for the input to return low.

The output of a channel in counter mode is always high.

To try the counter/timer mode:

One good way to see the counter in action is to have it compare the outputs of two clock channels. The counter will measure the difference between the successive pulses of the two channels.

1. On the General Settings, select Full (high resolution) clock mode.
Channels B and C will be our clocks. We’ll use odd delay and duration values for the clocks so that the difference between pulses will be obvious.
2. On Channel B’s properties tab, select Clock mode. Enter a delay of 54 ms and a duration of 79 ms. Click “Apply.”
3. On Channel C’s properties tab, select Clock mode. Enter a delay of 87 ms and duration of 85 ms. Click “Apply.”
4. On Channel A’s Logic tab, enter: $\text{inA} = \text{outB} \text{ and } \text{outC}$; (You’ll learn more about channel logic in the *Logic* section above).
5. On Channel A’s Properties tab, select Counter mode, with milliseconds as the units. Click “Apply.”

6. Lastly, click “Send Logic” on the Main screen to update the new logic for Channel A.

The time difference measurements should now appear, in rapid succession, in the “Last Measured Period” box on Channel A.



NOTE: The LC880 may be measuring time differences faster than it can display them all on the monitor.

Now, check “Hold measurement.” Note that the displayed measurement does not change. Un-check “Hold Measurement”, then check “Hold measurement” again to freeze another measurement.

Dynamically-Delayed Pulse Mode


The dynamic delay is the LC880’s *piece de resistance*. First, it measures the time difference between two pulses, like a counter. Next, it scales that difference by a user-defined factor. It subtracts an offset value, waits the duration of the dynamically calculated delay, then sends a pulse to the channel output.



NOTE: The dynamic delay is particularly useful for ballistics or gas dynamics experiments, in which the velocity of a projectile or wave must be calculated and used to time the triggering of an experimental device.

Dynamically Delayed Pulse mode includes several properties tabs.


Scaling



In the calculation of the delay, the measured delay between two pulses is scaled by a factor of 2^N , where $-16 < N < 16$, or $N/16$, where $0 < N < 256$. The former scaling is faster and a scaling of 1:1 is fastest. Scaling can typically be done at rates > 50 kHz.

The maximum delay between input pulses is 52 ms in hi-res mode (100 MHz) and 3.3 s in low-res mode (100/64). The maximum delay resulting from the scaling and offset calculations is respectively 21.47 s and 1374 s in hi-res and low-res modes. The delay must be positive.

Delay



Added Delay

Set an offset value to subtract from the scaled delay length. This is useful when accounting for camera shutters or other lag times in experimental equipment. The offset value can be positive (shortens the delay time) or negative (lengthens the delay time) and the resolution is set by the time base (e.g., 10 ns).

Duration of pulse

The duration of the output signal

Until Reset

The output will remain high after a trigger until the channel is reset

Inverted Output

Check to set the output high until a valid trigger. The duration of the pulse and other properties are the same as for the delayed trigger.

Triggering

The triggering options are the same as those for Delayed Pulse After Triggering mode:

Normal	The channel can retrigger an unlimited number of times
Skip N triggers	The channel triggers only after the specified number of input pulses (between 0 and 2,147,483,647).
Retrigger N times	Limits the number of retriggers allowed between resets (between 0 and 2,147,483,647). To allow a channel to be triggered again after reaching the limit, press the “Reset” button on the LC880 or choose Program>Reset Channels in Trigger.
Trigger on rising input	Sets the channel to be edge triggered. In other words, an input signal that stays high triggers the channel only once.
Trigger on input high	This option is not used in this mode.
Alert host computer following trigger	Sends a “ping” to the PC upon completion of the trigger cycle.

The front panel LED for a dynamically-delayed trigger will remain green until a valid input signal is detected. When an input signal is detected, the LED will glow red until the output signal is sent. At this point the LED will return to green if retriggering is allowed, or it will glow red to indicate that no more input signals will be accepted.

To try the dynamically-delayed pulse mode:

You can simulate the experiment above, using the trigger button to represent Sensor 1, and a delayed trigger in Channel B to represent the signal from Sensor 2. Connect your oscilloscope to the Output C to view the results.


1. On the General Settings screen, select 100/64 (low-resolution) mode for Channels C and D (and optionally for Channels A and B).
2. On Channel A’s Properties tab, select “Output=Input” mode.
3. On Channel A’s Logic tab, enter: inA = not in1;
4. On Channel B’s Properties window, select “Delayed Trigger” mode, with a delay and duration of 100 ms each. Setup for edge triggering and allow retriggering.
5. On Channel B’s Logic tab, enter: inB = outA; This logic will cause the delayed trigger to pulse 100 ms after a signal from Channel A.

6. On Channel C's Properties window, select "Dynamically Delayed Pulse" mode. Set up for scaling by 2^N , where $N = 1$, and set the Added Delay to $-200 \mu\text{s}$ and the duration to 100 ms. Allow unlimited retriggering.
7. Since you've made changes to multiple channels, go to the Main screen and click the "Set All Channels" button. All of your changes will now be stored in the LC880.

Now the test is ready. Depress and hold the trigger button. Output A will transition high, causing the delayed trigger in Channel B to begin its count. When Channel B sends its output pulse, the dynamic delay in Channel C will calculate the time difference (100 ms), multiply it by $32/16$, and subtract the $200 \mu\text{s}$ offset. 199.8 ms after receiving the signal from Channel B, Channel C will send a 100 ms output pulse.

Toggled Output Mode

This mode toggles its output state upon every input pulse. It functions like a flip-flop and is useful for a variety of logical and control purposes. The minimum input-pulse state duration is 10 ns and the edge-to-toggle delay is < 20 ns.



The only user settings for this mode are the initial state, which may be either low or high.

The front panel LED for the channel will glow green when the output is low or red when the output is high.


Dynamic Delay Compensation Mode

Many electronically triggered devices exhibit significant drift due to temperature changes, time, etc. This mode is designed to dynamically adjust trigger timing to compensate for drift in a laser or other electronics.

To use the mode, you'll enter the desired delay between the triggering of an event and the actual occurrence of the event (say, the light emission from a laser). Two signals are then mapped logically to the channel input: the trigger signal, and some indicator of the actual occurrence of the event (say, a photodiode triggered by the laser light). The channel then measures the actual delay between the trigger and response. It then adjusts the next trigger pulse accordingly to maintain the desired delay to within approximately 10ns. To reduce jitter from shot-to-shot timing fluctuations or noise you can determine the delay compensation by averaging over a series of trigger pulses rather than adjusting for each individual trigger.

Dynamically Delay Compensation mode includes several properties tabs.

Delay



Desired delay	The delay between the trigger and response signals that Trigger will attempt to maintain.
Duration of pulse	The duration of the output signal
Inverted Output	Check to set the output high until the next valid trigger.

Averaging

Delay Averaging	Determines how the delay compensation will be calculated. Choose None to calculate a new delay compensation time with each trigger pulse; or, choose to use an average compensation based on the last 2 to 256 trigger pulses. Averaging is performed with an infinite impulse response filter (i.e., a conventional one-pole low pass filter). The number of triggers sets the filter response time.
------------------------	--

Sweep

The delay between the trigger and response can be “swept” or incremented by a specified value each time the channel is triggered. The sweep value must be positive. The increment resolution is 5 ns. However, the actual delay is rounded to the nearest 10 ns at full internal clock. For example if the increment is set to 5 ns, the actual delay is incremented by 10 ns every second trigger event.

To try the dynamically-delayed compensation mode:

In this example an external delay generator is used to simulate a laser with the Trigger In as the laser trigger and the Delay Out as a signal from a photodiode monitoring the laser pulse. The delay within the delay generator will changed from 0 to 10 μ s, yet, regardless of this fluctuation, the laser will continue to fire 25 μ s after its trigger pulse.

Setup the LC880 as shown below:

Channel A: Free-running clock: High duration: 10.000 us; low duration: 990.000 us.

Channel B: Dynamic delay compensation: Adjust delay for 25.000 us between inputs.

Pulse output for 8.000 us.

Average 8 delays.

Trigger Input Logic: $\text{inB} = \text{outA} \text{ or } \text{in1}$;

Channel C: Delayed pulse: Delay after trigger 10.000 us then pulse output 8.000 us.

Triggering options: Trigger on rising edge. Unlimited retriggering.

Trigger Input Logic: $\text{inC} = \text{outB}$;

Channel D: Not used

Channel E: Delayed pulse: Delay after trigger 16.000 us then pulse output 10.000 us.

Triggering options: Trigger on rising edge. Unlimited retriggering.

Trigger Input Logic: $\text{inE} = \text{outA}$;

Channel F: Not used

Channel G: Not used

Channel H: Not used

Connect Channel E Out to the scope External Trigger.

Connect Channel C Out to the scope Channel 1 and Trigger In of delay generator.

Connect Delay Out of delay generator to scope Channel 2 and in1 of LC880.

Set up the delay generator so that it has a variable delay between 0 and 10 μs to simulate the drifting delay between when the laser is triggered and when it actually fires. A few considerations must be noted:

1. The pulse out from B cannot overlap the feedback pulse. Since pulse widths in the LC880 must be $\geq 7.695 \mu\text{s}$, a second delay unit must be placed between the LC880 and the laser if the time between triggering and firing of the laser is $< 7.695 \mu\text{s}$. Channel C serves this purpose in this example.
2. The minimum delay for Channel C should exceed the pulse width used in Channel B.
3. For the feedback circuit to function, the dynamic delay compensation time must be greater than the sum of the Channel B pulse width ($\geq 7.695 \mu\text{s}$), the Channel C delay ($\geq 7.695 \mu\text{s}$) plus approximately 7 μs to perform the calculation, which is approximately 23 μs .

RUNNING EXPERIMENTS

Once you've created a program you're ready to set up and perform your experiment.

Storing and Recalling Experiment Files

As mentioned earlier, an experiment file is a document which resides on your PC, containing all of the logic, properties, and settings profiles for a particular experiment. Experiment files bear the .trg extension.

To save an experiment, choose "Save As" from the File menu. Name the file, then store it on your hard drive in an accessible place.

To reuse an experiment file, choose File>Open, and locate the file. When you click "OK," the file's settings will appear in Trigger. Click "Send Properties," then click "Send Logic" to download the settings to the LC880. You can also choose "Program All Channels" and "Program Logic" from the Program menu.

To print a text summary of the settings in the current experiment, choose File>Print.



NOTE: Good file naming protocol will help your experiments go more smoothly. For example, if you are saving multiple versions of a combustion experiment, store all files in the same directory, with the names "combustion1.trg", "combustion2.trg", etc. Or, better yet, use descriptive names such as "combust - 50us delay.trg", "combust - 60us delay.trg", etc.

Front Panel Controls and Indicators

The front panel of the LC880 is designed to help you run experiments without computer control. This functionality allows you to program the unit at a workstation, then to move the unit to the experimental environment.

Channel LEDs

As mentioned earlier, the LEDs for each channel provide information on each channel's status:

1. A green light on a channel typically indicates that the channel is waiting for an input trigger.

When a channel is in Fixed Output, Clock Pulse, Output=Input or Toggled timing modes, the green light indicates that the output is currently LOW for that channel.

2. A red light typically indicates that the channel is in the process of triggering.

When a channel is in Fixed Output, Clock Pulse, Output=Input or Toggled timing modes, the red light indicates that the output is currently HIGH for that channel.

3. When a channel indicator is steadily off, the channel is disabled (see *Enable/Disable Channels* below). The channel indicator also turns off during delays and while waiting for the input to return low for an edge trigger.

The following chart provides an overview of the front panel channel LEDs:

Mode	Off	Green	Red
Fixed Output	--	Output is Low	Output is High
Output = Input (Passive)	--	Output is Low	Output is High
Clocked Pulse Stream	--	Output is Low	Output is High
Delayed Pulse After Trigger	In delay period	Waiting for trigger; if retriggering is on, light will glow green once input goes low	In duration period; if retriggering is off, light remains red after duration period
Validate Trigger (Noise Suppression)	--	Waiting for trigger; if retriggering is on, light will glow green once input goes low	In duration period; if retriggering is off, light remains red after duration period
Counter	Count data is being reported, waiting for input to go low	Waiting for first trigger pulse	Waiting for second trigger pulse
Dynamic Delayed Pulse	--	Waiting for trigger; if retriggering is on, light will glow green once input goes low	In duration period; if retriggering is off, light remains red after duration period
Toggle	--	Output is Low	Output is High



NOTE: Channel LEDs provide important information about the channels, but an oscilloscope remains a useful tool for diagnosing channel outputs, particularly for high-speed experiments.

ALL Channel LED

The All Channel LED shows the overall status of the channels.

A green LED indicates that the external inputs and outputs are enabled and functioning normally

A red LED indicates that the external inputs and outputs are frozen. (See *Running Experiments: Disable/ Enable channels* below).

Selecting a Channel from the Front Panel

Each front panel function requires that you select a channel (or ALL channels) and a function.

To select a channel, click the channel button. The channel LED will glow brightly, as will the available Function buttons. The ALL channels button will function similarly to select all channels. You can now perform one of the front panel functions on that channel.

In the image below, the F channel button has just been pressed. The buttons for the four available functions (Enable, Disable, Trigger and Reset) glow yellow.



Performing a Front Panel Function

Front panel functions can be accessed in one of two ways:

- select a channel (or ALL channels), then choose a function
- select a function, then choose a channel (or ALL channels).

There is no significant difference between these two approaches.

In the image below, the Enable button has just been pushed. The channel buttons and the ALL channel button are lit. Pressing one of these buttons will now enable the channel(s).



Enable/Disable ALL

Occasionally you'll need to abort an experiment in progress, due to equipment problems, etc. The quickest way to shut down an experiment is to disable the inputs and outputs of the LC880. To do so, click the "Disable" button, then click a particular channel button, or the ALL channels button. The inputs and outputs for all channels will be immediately frozen.

Following a Disable command, you can re-enable the inputs and outputs by clicking the "Enable" button, then clicking a channel button or the ALL channels button.

Enable/Disable Channel

It is sometimes convenient to disable or suspend a channel temporarily. The Enable/Disable channel functions provide this capability from the front panel.

You can also disable a channel, or all channels, from within Trigger (see *Trigger Software Interface: General Settings: Disable/Enable Channels*).

Trigger a Channel

This function lets you simulate an input trigger on a channel, regardless of the channel's logic settings. This is a useful way to initiate a timing sequence without actually performing the live experiment. You can also trigger a channel from within the Trigger software (See *Channel Properties: Force Trigger* above).



NOTE: You cannot trigger ALL channels, only individual channels.

Reset

The Reset button returns the selected channel or ALL channels to the last programmed settings. It resets and reinitializes retriggering counters and the delays of any channels that are undergoing a delay sweep.

Store/Recall Settings Profiles

As you know, you can store 64 different “Settings Profiles” in the LC880. A Settings profile includes properties, variables, and logic for each of the eight channels. Storing Settings profiles lets you switch between settings while running an experiment, to test different timing criteria, initial conditions, etc., without re-plugging or reprogramming.

Settings profiles are enumerated AA, AB..., BA, BB... HH, and they are accessible under those letters from the LC880.

To store the current settings in a Settings profiles:

1. Click the Store button
2. Click the two-letter combination for the Settings profile.


Similarly, to recall a stored Settings profile:

1. Click the Recall button
2. Click the two-letter combination for the Settings profile.

You can recover these settings even after the unit has been turned off or disconnected. The most recently used setting is automatically loaded when the LC880 is powered up. If you attempt to recall an invalid setting, the power LED flashes red and the system retains its previous setting.

To recall a setting profile from within Trigger:

1. Click the “Recall Setting...” from the Program menu.
2. Click on or type in the desired setting (case insensitive).
3. Click “Recall.”



To store a setting profile from within Trigger:

1. Click the “Store Setting...” from the Program menu.
2. Click on or type in the desired setting (case insensitive).
3. Click “Store.”

GLOSSARY

Apply	Sends a channel's current timing settings to the LC880.
Baud Rate	The speed of communication via the Com port.
Channel	One of eight internal signal paths in the LC880.
Channel Input	The signal which feeds one of the eight channels. A channel input can be triggered by signals from any combination of signals from the input jacks and/or channel outputs.
Channel Settings	The logic and timing settings which determine how a channel functions.
Com Port	A serial port on a PC computer. The LC880 communicates with the computer via an RS-232 cable connected to a Com port.
Controller	The LC880 hardware unit.
Counter/Timer Mode	This timing mode determines the time delay between two separate input signals.
Delay	The time between the arrival of an input trigger and the release of an output pulse.
Delayed Pulse After Trigger	In this timing mode, a channel waits for an input signal, waits a specified delay time then sends a pulse to the channel output.
Duration	The length of an output pulse.
Dynamic Delay Compensation	This mode is designed to dynamically adjust trigger timing to compensate for drift in a laser or other electronics.
Dynamically Delayed Pulse	This timing mode measures the time difference between two pulses, scales that difference by a user-defined factor, waits the duration of that calculated delay, then sends a pulse to the channel output.
Edge Triggering	This channel setting ensures that an input signal which stays high will only trigger the channel once.
Enable/Disable I/O	A function which disables or suspends a channel (or all channels) temporarily. The function is available on the LC880 front panel and on the channel settings screens.

External Clock	The “Ext Clk” input on the back of the LC880 allows the time base to be controlled by an external source.
Firmware	The internal code of the LC880 controller. The firmware, as well as the Trigger control software, can be field-upgraded. Latest upgrades are available at www.labsmith.com .
General Settings	The settings which govern the function of <i>all</i> channels.
Global Variable	A shorthand symbol for a logical assignment that may be used repeatedly throughout a program.
High Resolution	A timing mode in which the clock operates at 100 MHz.
Implement Logic	A menu function which sends the logic for all channels to the LC880.
Input	One of eight inputs on the back of the LC880. Signals from these inputs can be combined to trigger any channel.
Local Variable	A preset variable that is used in the logical assignments that control channel triggering.
Lock/Unlock Panel	A function that disables all front panel functions (including the power button) to avoid tampering/accidental changes when the LC880 is embedded in an experiment.
Logic	The set of assignments that determines when a channel will trigger.
Low Resolution	A timing mode in which the clock operates at 100/64 MHz, or 1.56 MHz.
Noise Suppression	Also called “Validated Trigger Input” mode. This timing mode eliminates false triggers due to noise on the input signal by requiring the signal to remain high (or low) for a specified duration.
Output	One of eight physical output jacks on the back of the LC880. The output signals are determined by the eight timing channels.
Parsing	Parsing occurs when you click the “Apply” button. The Trigger software checks your logic assignments and changes them to code for the LC880 controller.
Predefines	A set of standard variables and operators that can be used in logic assignments.

Program All Channels	A function which sends all logic and timing settings to the LC880.
Programmable Logic Controller	A device which combines logic and timing processors for precise control of physical experiments.
Recall	Re-loads a stored Settings File.
Reset	Returns all channels of the LC880 to their waiting condition. Retriggering Counters and delay settings are reset as well.
Retrigger	A channel setting that allows the channel to be triggered more than once, if subsequent trigger pulses are received.
Scaling	In the Dynamically Delayed Trigger timing mode, scaling determines the ratio of the distance between the two input sensors to the distance of the target.
Settings File	A PC file (with a .trg extension) that contains all the logic and timing settings for an experiment. Up to 256 Settings Files can be stored and recalled from the LC880.
Skip Triggers	A setting that allows a channel to trigger only after the specified number of input pulses.
Store	Load a Settings File into the memory of the LC880. The settings can be recalled later.
Sweep	In the Delayed Pulse After Trigger timing mode, this settings facilitates time-delay scanning experiments in which the delay of the pulse relative to the trigger must be incremented by a specified value each time the channel is triggered.
Time Base Settings	These General Settings determine the timing resolution of the LC880. Higher resolution allows 10ns timing resolution, while low resolution allows delays over 1300 seconds.
Timer	Counter/Timer mode determines the time delay between two separate input signals and reports it back to the PC.
Timing Mode	The eight functions that each channel can perform.
Toggle Mode	This mode toggles its output state upon every input pulse.
Trigger	An input signal that initiates a channel function.
Trigger (software)	The control software for the LC880.

Trigger Validation Mode

Also known as “Noise Suppression” mode. This timing mode eliminates false triggers due to noise on the input signal by requiring the signal to remain high (or low) for a specified duration.

Update Firmware

This menu function allows you to upgrade your LC880’s firmware to the latest version, available at www.labsmith.com.

Trigger and LC880 are trademarks of LabSmith.

National Instruments is a registered trademark of National Instruments, Inc.

LabView is a trademark of National Instruments, Inc.

Microsoft, Windows 2000 and XP are registered trademarks of Microsoft Corporation.

Adobe and Acrobat are registered trademarks of Adobe Systems, Inc.

Trigger software ©1996--2008 LabSmith.

This manual ©2008 LabSmith. No part of this document may be reproduced or distributed without the consent of LabSmith.