


Reynolds Number Worksheet

The **Reynolds Number (Re)** is a dimensionless number that can be used to predict flow patterns in different fluid flow situations, such as air over an aircraft wing or liquid in a pipe. Reynolds number calculations are often used to design products or experiments. For instance, if you were testing the efficiency of a small-scale model of a car in a wind tunnel you would design your experiment so the full scale and model have similar Reynolds Numbers.


For flow in a tube or pipe, the Reynolds Number can be used to predict whether the flow will be **laminar** or **turbulent**.

Definitions
dimensionless number: a number representing a physical property that has no unit scale
Reynolds Number: the ratio of inertial forces to viscous forces
turbulent flow: flow in which the fluid undergoes irregular fluctuations
laminar flow: flow in which the fluid moves in smooth paths or layers
Stokes Flow: fluid flow where inertial forces are small compared with viscous forces

Re < 2000




Re > 4000




In microfluidics, the Reynolds Number is typically very low, often much less than 1, and this is a kind of extreme laminar flow known as **Stokes Flow**. In Stokes Flow, fluids will stay laminar even when flowing around corners or obstacles. This can affect things like heat transfer and reaction rates, as flow through a microfluidic channel typically only mixes through diffusion.

Re ≈ 0.01



Re ≈ 20



Re ≈ 2,000



The Reynolds number depends on the fluid density, velocity and viscosity as well as the dimensions of the channel. Since the Reynolds Number is a **dimensionless number**, all units (m, kg, etc.) will cancel out for the final value.

$$Re = \frac{\rho v d}{\mu}$$

ρ = density

v = velocity

d = diameter of channel

μ = dynamic viscosity

For water at ~20 C

$$\rho_{water} = 1000 \text{ kg/m}^3$$

$$\mu_{water} = 1.0 \text{ N-s/m}^2 = 1.0 \text{ kg/m-s} \quad (1 \text{ N} = 1 \text{ kg-m/s}^2)$$

Example: Will water flowing through a 1 cm diameter circular pipe at 10 cm³/s be laminar or turbulent?

What we know:

density of water: $\rho_{water} = 1000 \text{ kg/m}^3$

dynamic viscosity of water: $\mu_{water} = 1 \frac{\text{kg}}{\text{m}\cdot\text{s}}$

flow rate in the channel: $F = 10 \text{ cm}^3/\text{s}$

channel diameter: $d = 1 \text{ cm}$

$$Re = \frac{\rho v d}{\mu}$$

You have all needed values to plug into the Reynolds Number equation except the velocity (v). You know the flow rate (volume/time), but you need to know the velocity (distance/time).

Step 1: Find the velocity in m/s:

$$v = \frac{F}{A}$$

Where F = flow rate and A = area.

$$A_{circle} = \pi r^2 = \pi \left(\frac{d}{2}\right)^2 = \pi \left(\frac{1 \text{ cm}}{2}\right)^2 = 0.78 \text{ cm}^2$$

$$v = \frac{F}{A} = \frac{10 \frac{\text{cm}^3}{\text{s}}}{0.78 \text{ cm}^2} = 12.8 \frac{\text{cm}}{\text{s}} \times \frac{1 \text{ m}}{100 \text{ cm}} = \mathbf{0.128 \frac{m}{s}}$$

What's the difference between flow rate and velocity?

flow rate is the *volume* of fluid that moves past a surface per unit time

velocity (or speed) is the *distance* the fluid moves per time.

To calculate the velocity from the flow rate, use the equation:

$$velocity = \frac{\text{Flow rate}}{\text{Cross sectional area}}$$

Step 2: Convert the diameter into meters:

$$d = 1 \text{ cm} \times \frac{1 \text{ m}}{100 \text{ cm}} = 0.01 \text{ m}$$

Step 3: Plug into the Reynolds Number equation:

$$Re = \frac{\rho v d}{\mu} = \frac{(1000 \frac{\text{kg}}{\text{m}^3})(0.128 \frac{\text{m}}{\text{s}})(0.01\text{m})}{(1.0 \frac{\text{kg}}{\text{m}\cdot\text{s}})}$$

Step 4: Make sure all of your dimensions cancel out and solve for Re.

$$Re = 1.28 \frac{\left(\frac{\text{kg}}{\text{m}^3}\right)\left(\frac{\text{m}}{\text{s}}\right)(\text{m})}{\left(\frac{\text{kg}}{\text{m}\cdot\text{s}}\right)}$$

Re = 1.28

This flow will be laminar.

Problem 1: You have water flowing in a 2.5 mm diameter circular channel at a flow rate of 10 $\mu\text{l/s}$.

a. What is the fluid velocity in m/s?

Hint: 1 $\mu\text{l} = 1 \text{ mm}^3$

b. What is the Reynolds Number?

$$Re = \frac{\rho v d}{\mu}$$


Remember:

$$\rho_{\text{water}} = 1000 \text{ kg/m}^3$$

$$\mu_{\text{water}} = 1.0 \text{ N-s/m}^2 = 1.0 \text{ kg/m-s}$$

The Reynolds Number equation assumes a circular cross section in the channel. You can calculate the *hydraulic diameter* of a non-circular chip with the equation:

$$d_H = \frac{4A}{P}$$
$$= \frac{4 \times (\text{Cross sectional area})}{\text{Perimeter}}$$



Problem 2: The main channel in the Y-microfluidic chip included in the Basic Microfluidic Education Kit has a rectangular cross-section with a width of 0.4 mm and depth of 0.1 mm.

- What is the channel hydraulic diameter?
- If the flow rate across the chip is 100 $\mu\text{l}/\text{min}$, what is the velocity in meters per second?

Remember: 1 $\mu\text{l} = 1 \text{ mm}^3$

- Calculate the Reynolds Number for the channel, assuming water flow. Don't forget, the Reynolds Number is a dimensionless number, so all units will cancel.
- With this Reynolds Number, would you expect the fluids from the two channels to mix? Explain.

Reynolds Number Worksheet - Solutions

Problem 1: You have water flowing in a 2.5 mm diameter circular channel at a flow rate of 10 $\mu\text{l/s}$.

a: What is the fluid velocity in m/s?

$$v = \frac{F}{A}$$

$$F = \text{volume flow rate} = 10 \mu\text{l/s} = 10 \text{ mm}^3/\text{s}$$

$$A = \text{area} = \pi r^2 = \pi \left(\frac{d}{2}\right)^2 = \pi \left(\frac{2.5}{2}\right)^2 = 4.9 \text{ mm}^2$$

$$v = \frac{F}{A} = \frac{10 \frac{\text{mm}^3}{\text{s}}}{4.9 \text{ mm}^2} = 2.04 \frac{\text{mm}}{\text{s}} \times \frac{1 \text{ m}}{1000 \text{ mm}} = \mathbf{0.002 \frac{m}{s}}$$

b. What is the Reynolds Number?

$$\text{Re} = \frac{\rho v d}{\mu} = \frac{(1000 \frac{\text{kg}}{\text{m}^3})(0.002 \frac{\text{m}}{\text{s}})(2.5 \text{ mm})}{(1.0 \frac{\text{kg}}{\text{m s}})}$$

$$\text{Re} = \frac{(1000 \frac{\text{kg}}{\text{m}^3})(0.002 \frac{\text{m}}{\text{s}})(2.5 \text{ mm})(\frac{1 \text{ m}}{1000 \text{ mm}})}{(1.0 \frac{\text{kg}}{\text{m s}})}$$

$$\mathbf{\text{Re} = 0.005}$$

Problem 2: The main channel in the Y-microfluidic chip included in the Basic Microfluidic Education Kit has a rectangular cross-section with a width of 0.4 mm and depth of 0.1 mm.

a: What is its hydraulic diameter?

$$d_H = \frac{4A}{P}$$

$$A = \text{cross sectional area} = w \times D = 0.4 \text{ mm} \times 0.1 \text{ mm} = 0.04 \text{ mm}^2$$

$$P = \text{perimeter} = 2 \times (w + D) = 2(0.4 \text{ mm} + 0.1 \text{ mm}) = 1 \text{ mm}$$

$$d_H = \frac{4A}{P} = \frac{4 \times 0.04 \text{ mm}^2}{1 \text{ mm}} = \mathbf{0.16 \text{ mm}}$$

b: If the flow rate on the chip is 100 $\mu\text{l}/\text{min}$, what is the velocity in meters per second?

$$v = \frac{F}{A}$$

$$F = \text{volume flow rate} = 100 \mu\text{l}/\text{min} = 100 \text{ mm}^3/\text{min}$$

$$A = \text{cross sectional area} = 0.04 \text{ mm}^2$$

$$v = \frac{F}{A} = \frac{100 \frac{\text{mm}^3}{\text{min}}}{0.04 \text{ mm}^2} = 2500 \frac{\text{mm}}{\text{min}} \times \frac{1 \text{ min}}{60 \text{ s}} \times \frac{1 \text{ m}}{1000 \text{ mm}} = \mathbf{0.041 \frac{\text{m}}{\text{s}}}$$

c: Calculate the Reynolds Number for the channel, assuming water flow.

$$\text{Re} = \frac{\rho v d_H}{\mu} = \frac{\left(1000 \frac{\text{kg}}{\text{m}^3}\right) \left(0.041 \frac{\text{m}}{\text{s}}\right) (0.16 \text{ mm})}{\left(1.0 \frac{\text{kg}}{\text{m s}}\right)}$$

$$\text{Re} = \frac{\left(1000 \frac{\text{kg}}{\text{m}^3}\right) \left(0.0041 \frac{\text{m}}{\text{s}}\right) (0.16 \text{ mm}) \left(\frac{1 \text{ m}}{1000 \text{ mm}}\right)}{\left(1.0 \frac{\text{kg}}{\text{m s}}\right)}$$

$$\mathbf{\text{Re} = 0.00066}$$