

4AM01 Sensor Manifold

- ▶ Simultaneous monitoring and control of up to four uProcess™ sensors
- ▶ Interfaces with uProcess software for coordinated control




Analog sensor manifold controls up to 4 uPS pressure sensors and/or uTS temperature sensors using uProcess software.

Specifications

Interface to sensors	6 conductor (1 mm pitch) same-side flat flex ZIF cable. 4AM01 includes 4 flat flex ribbon cables: 1 x 8", 2 x 10", 1 x 12"
EIB200 Automation Interface	Required for control of 4AM01 via uProcess™ software; includes software, LabVIEW™ drivers, and SDK
Interface to EIB200 or breadboard	2.00 mm pitch, 20 pin, dual row receptacle
Case material	Plated aluminum
Dimensions	2.5" x 0.7" x 1.1" (64 mm x 18 mm x 28 mm)

Dimensions



3D model available upon request