

4VM02 Valve Manifold

- ▶ Simultaneous control of up to four uProcess™ automated valves
- ▶ Interfaces with uProcess software for coordinated control




The 4VM02 manifold connects and controls up to four LabSmith uProcess™ automated valves (AV201, AV202, AV303, or AV801). Requires EIB200 automation interface for use with uProcess software. Mounts to uProcess breadboards (uPB-5, uPB-8, uPB-12, or iBB).

Specifications

Interface to valves	6 conductor (1 mm pitch) same-side flat flex ZIF cable. 4VM02 includes 4 flat flex ribbon cables: 1 x 8", 2 x 10", 1 x 12"
EIB200 Automation Interface	Required for control of 4VM02 via uProcess™ software; includes software, LabVIEW™ drivers, and SDK
Interface to EIB or breadboard	2.00 mm pitch, 20 pin, dual row receptacle
Case material	Plated aluminum
Dimensions	2.3" x 0.4" x 1.1" (58 mm x 11 mm x 28 mm)

Dimensions



3D model available upon request