


EIB200 Automation Interface

- ► For monitoring and controlling uProcess[™] devices
- ► Control over 100 devices
- USB computer interface
- ▶ uProcess[™] software, LabVIEW[™] drivers, SDK included




Specifications	EIB200
Computer interface cable (included)	Micro USB 2.0
Case material	Conductive ABS
Power supply (included)	12 V/ 2.0 A
Power jack connector	2.1 mm x 5.5 mm
Control and data acquisition	uProcess ™ Software, LabVIEW™ drivers, SDK (C++)
Interface to uProcess devices and breadboard	2.00 mm pitch, 20 pin, dual row receptacle

Accessories	Description
uProcess breadboard uPB-05, uPB-08, uPB-12, or iBB	Optional platform for mounting EIB200 and other uProcess and CapTite microfluidic components



EIB uProcess Interface

Specifications



26cm

3D model available upon request